

# SIEMENS

## Uroskop ACCESS

SP

### Update Instructions

SP015/03/S

Title: New arm holder for table

Reason for Update:

☒ Safety

Urgency:

☒ Immediate

☐ Within

Update material required?

☒ Yes

☐ No

Materials free of charge?

☒ Yes

☐ No

Return of parts?

☒ Yes

☐ No

Estimated completion time: 0,5 h

Number of CSE's:1

Customer application training? ☐ Yes

☒ No

#### Systems/Products affected/System identifying IVK

Name	Material No.	Serial No.
Uroskop ACCESS LA	57 56 122	1057 and 1060 - 1070
Uroskop ACCESS RA	57 56 130	2025 - 2032

Remark: n.a.

#### Components affected/to be modified

Name	Material No.	Serial No.	Component status affected
n.a.	n.a.	n.a.	n.a.

Remark: n.a.

Chg. Ref. No.:

127221

Name:

Kölnner

Dept.:

CS PS 24

© Siemens AG 2003

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

Print No.: SPL5-330.896.06.01.02

Doc. Gen. Date: 11.03

Replaces: n.a.

Version 2.1; January 14, 2003

Page 1 of 6

**Document Revision Level**

This document corresponds to the version/revision level effective at the time of system delivery. Revisions to hardcopy documentation are not automatically distributed.

Please contact your local Siemens office to order current revision levels.

**Disclaimer**

The installation and service of equipment described herein is to be performed by qualified personnel who are employed by Siemens or one of its affiliates or who are otherwise authorized by Siemens or one of its affiliates to provide such service.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens or one of its affiliates are directed to contact one of the local offices of Siemens or one of its affiliates before attempting installation or service procedures.

**Systems/Products Affected**

Uroskop ACCESS with System no.:

57 56 122 Uroskop ACCESS LA, serial no. 1057 and 1060 - 1070

57 56 130 Uroskop ACCESS RA, serial no. 2025 - 2032

**Reason for the Update**

The arm shield has to be changed because of the potential risk of injury.

**Prerequisites**

n.a.

**Special Tools / Documents**

n.a.

## Ordering Information

The following update kit has to be ordered from CSML (SAP Distribution Channel, factory 2050) or BU Logistics (factory BU) (only with an order via BU Logistics):

### Update Kit 77 57 482

List the system Serial Number on the order!

## Contents of the Update Kit

Update kit **77 57 482** contains the following parts:

Pos.	Quan.	Material No.	Name
1	1	73 09 797	Arm Shield
2	1	SPL5-330.896.06.01.02	Update Instructions
3	1	SPL5-330.623.06.01.01	Supplement to Instruction
4	1	SPL5-330.623.06.01.02	Supplement to Instruction
5	1	SPL5-330.623.06.01.11	Supplement to Instruction
6	1	SPL5-330.623.06.01.12	Supplement to Instruction
7	1	SPL5-330.623.06.01.24	Supplement to Instruction

Tab. 1

## Return of Parts

n.a.

**Work Steps**

- Attach the arm shield as shown in the "Operating Instructions Urooskop ACCESS", Supplement - Arm Shield, included with delivery.
- Inform the staff about how to handle the arm shield.

**Final Check**

- Function check of the system according to the modifications that have been made.

**Customer Information**

The customer information SPL5-330.623.06.01.xx is included in delivery.

**Final Work Steps**

- Update the system documentation.  
Update the operating instructions.  
Fill out, and if needed, make a copy of the attached "Completion Protocol/ Update Completion Form" and file it in the corresponding System Binder/User Handbook.
- Updates that have already been completed prior to publication of this UI must also be reported.

**Changes to Previous Version**

n.a.

## Completion Protocol

The update with the number **SP015/03/S** has been completed.

Material number: .....

Serial number: .....

Customer: ..... Functional Location: .....

Customer No.: .....

Name ( CSE ): ..... Telephone: .....

Country: ..... Location: .....

Date: ..... Signature: .....

Remark: .....

.....

<b>NOTE</b>
-------------

**After completing the update, make a copy of this page, fill it out and file it in the corresponding System Binder/User Handbook.**

This page intentionally left blank.